



# Insulin Travel Pen Wallet

Keep your insulin cool  
while on the go



## FACT BOX

**Product Code:**  
MD117

**Product Category:**  
Home Health

### Contents:

Wallet that can hold x2 insulin pens  
(pens not included)



No need for  
refrigeration



Ideal for people  
on the go



Activates in  
40 seconds

Medicare® Insulin Pen Travel Wallet is specially designed to offer diabetes sufferers and other medication users the ultimate accessory for the safe and convenient storage/ transportation of their insulin pens, vials and other diabetic supplies they might need at correct and safe temperature without the need of refrigeration. It can be reactivated by placing in water when required.

### FEATURES:

- Patented fabric technology to keep insulin cool
- Can be re-activated multiple of times
- Great for people on the go
- Quick and simple to use

### DIRECTIONS:

1. Soak pack in cold water for 1-2 minutes or hold and turn vigorously under running tap water for 30-40 seconds.
2. Squeeze out any excess water and pat dry with a towel.
3. Insert items to be cooled and place pouch back inside wallet.
4. Repeat every 3 to 5 days to reactivate.

### CONTRAINDICATIONS AND PRECAUTIONS:

If not using the pouch, remove the pockets and allow to air dry.



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid